

# Alternative analysis

TECHNICAL PROFILE

## An internet-enabled analytical NIR system provides real-time results and supply chain integration without the need for specialist expertise

by Suzi Fraser Dominy and Emily Buckley

In the thirty years since it was first introduced, NIR — Near-Infrared Technology — has become a household word in the grain based industries: flour millers rely on this rapid and accurate method of analysis to determine the ash, protein and moisture content of flour; feed millers analyze raw materials for crude protein, fat and starch to optimize formulations and to assure the quality of finished feeds leaving the plant; while grain handlers depend on it to perform load-by-load protein checks, to set prices and make possible segregation of grains and oilseeds.

Yet despite its versatility, speed and accuracy, the technology is still not fully exploited, since it requires capital outlay, skilled operation and time and expense in initial and on-going calibration.


To overcome these drawbacks, an Internet-enabled NIR system has been developed by Cognis Corporation. Its unique QTA system, which stands for Quality Trait Analysis, takes a different approach to NIR systems as it supplies its clients with an FT-NIR machine and simply charges a fee per analysis, while offering the ease of remotely monitoring a client's unit and

### Grain cooling

## GRANIFRIGOR®

... the most natural way of grain preservation by cooling

- \* Protection against insects and microbes
- \* Without chemical treatment
- \* Short amortisation period
- \* Low energy costs
- \* Independent of ambient weather conditions



formerly Sulzer-Escher Wyss GmbH

**Axima Refrigeration GmbH**  
 Kemptener Str. 11-15 · 88131 Lindau, Germany  
 fon +49 8382 706-1 · fax +49 8382 706-410  
[www.axima.eu.com/refde](http://www.axima.eu.com/refde) · [info-refde@axima.eu.com](mailto:info-refde@axima.eu.com)

*Axima is a member of the eurammon initiative.*  
 For more information about eurammon see [www.eurammon.com](http://www.eurammon.com)

Circle No. 70 on the Reader Service Card or see page 72.

### Simple. Controllable. Selective. The UMS Discmill® type MHA600D



The Discmill® is simple to operate, simple to install and simple to maintain - as a result of its sophisticated design.

The Discmill grinds cereal grains controllably. A micrometer adjustment of the grinding gap allows it to produce wholewheat, rye flour, corn meal and many other cereal products to precise granulation specifications.

The Discmill grinds mill stocks selectively. It will produce low ash flour from a typical head reduction feed, just the way a roller mill does.

Consider the alternative from:

## UNITED MILLING SYSTEMS



HEAD OFFICE:  
 BANEMARKSVEJ 50 A • DK-2805 BRØNDBY • DENMARK  
 TEL. +45 70 22 27 11 • FAX +45 70 22 27 12 • E-MAIL: [umsmail@ums.dk](mailto:umsmail@ums.dk)  
 WEB: [www.ums.dk](http://www.ums.dk) • [www.umsamerica.com](http://www.umsamerica.com)

US OFFICE: TEL. +1 763 631 3003

UK OFFICE: TEL +44(0) 1995 672 556

Circle No. 71 on the Reader Service Card or see page 72.

June 2003 / **World Grain** / [www.World-Grain.com](http://www.World-Grain.com)



# Rice Processing Equipment.

From a highly reputed manufacturer.

Top quality rice comes from CHEN SAN FUNG rice processing machinery.

Vertical Rice Whitening Machine  
AMV-60 FTV-60 FTV-60

Rice Polisher  
MRP-40

Paddy Separator  
KSP-600

Paddy Husker  
CArBH-10M



**CHEN SAN FUNG MACHINERY INDUSTRIAL CO., LTD.**

No. 179, Sec. 3, Chung Shan Rd., Tantz, Taichung, TAIWAN.

Tel.: +886-4-2533-5135 Fax.: +886-4-2533-8855

E-mail: CSFM@ms8.hinet.net

http://www.CSFM.com.tw

Circle No. 74 on the Reader Service Card or see page 72.

A Fully Automated Revolutionary Rice Sorting Technology

**DEVELOPED BY YAMATEC**

**Unique in the World**

*Automatic operation with no manual interruption.*

Grain Art Sorter:  
With Two or Four  
Sensors

**Channel  
Options:**  
5ch,  
10ch,  
20ch,  
50ch,  
100ch,  
150ch

Quality Control Machine  
QCM-200



High Precision ACTIVE Electronic Sensor  
( Exceeds CCD Resolution )  
ART-150



**Agricultural Applications:**

White Rice  
Brown Rice  
Wheat  
Peanuts  
Sesame Seeds  
Soy Beans  
Red Beans  
Mung Beans  
Corn, etc.

**Industrial Applications:**

Gems  
Pearls  
Glass  
Plastics, etc

*Loop-back quality check during sorting process*

Manufacturer:

**Active KT Co., Japan,**

3-9-15 Minami-Chou, Chuou-Ku, Chiba City,  
Chiba, 260-0842, Japan



Worldwide Distribution and Service Center:

**Taiwan Yamakawa Technology Co., Ltd.**

3F, No. 23, Chung Shan N. Rd. Sec 2, Taipei, 104, Taiwan.

Tel.: +886-2-2561-6918 Fax: +886-2-2561-6977

Email: yamatec@ms47.hinet.net Website: www.yamatec.com.tw

Circle No. 75 on the Reader Service Card or see page 72.

maintaining all calibrations. Because Cognis monitors calibrations from off-site via an Internet connection, it is possible to ensure NIR units in separate locations have consistent results based on client-specified customizations.

**EASY TO USE**

At a client's facility, the QTA system comprises one or more NIR units, which is connected to a notebook computer that uses the Internet to access the Cognis-maintained central processor.

Calibrations are developed with clients for analyzing traits of materials that meet their specific needs, and calibration models are stored in the central processor.

Performing an analysis is simple, and users do not have to know how to choose NIR measurement parameters or how to develop/maintain calibration models. The user only has to go to QTA's website (www.QTA.com), select the test material from a drop-down list specific to each client

(such as canola seeds, wheat, soybean, sunflower, corn, etc.) and the traits to be analyzed (such as linolenic acid, moisture, protein, etc.). Users then follow step-by-step instructions to measure the sample's NIR spectrum. The obtained spectrum and the user input are sent to the central processor via the Internet connection.

The central processor uses the appropriate specified calibration models (according to the user input in earlier steps) and computes the sample for analysis. The results are then returned to the client.

Analysis data, along with sample identification and tracking information input by the customer, is stored in a secure central database, so that users can access the data to create reports or show records of sampling times. The stored information can be transferred, used and leveraged throughout an organization, providing ready access to comprehensive product quality information that could facilitate pricing and enhance marketing.



Cognis developed a Chingometrics centralized calibration system to increase consistency, eliminate the need for on-site recalibration and to provide simple operation.

**EQUIPMENT**

The FT-NIR sampling system is a rotating cup with a glass bottom and a diameter of about 100mm. The NIR beam illuminates an area approximately 20mm in diameter at one time.

When starting a measurement, the cup starts to rotate and a tumbling device also moves the sample material from the

# Seven Steps to Success

7 Modules covering a full spectrum of **Flour Milling Training...**

Hygiene Health and Safety

Wheat and the Screenroom

Mill Processes and Performance

Materials Handling Storage and Distribution

Flour

Power and Automation

Flour Milling Management

**Take control of your career and future**

For a full syllabus please contact nabim

21 Arlington Street  
London  
SW1A 1RN  
UK

Tel +44 207 4932521  
Fax +44 207 4936785  
www.nabim.org.uk  
e-mail: info@nabim.org.uk

Work towards the internationally recognised **City and Guilds Certificates in Flour Milling**

An indispensable overview of the industry for all employees

**Provides millers with an underpinning knowledge and understanding of the milling industry**

Studied every year by hundreds of millers worldwide...Europe, Africa, Asia, Australasia and North America

**Flexible training programme Supported by industry experts**



bottom layer to the upper layers. Therefore, this system can be used to analyze inhomogeneous solid materials including powder, grains, pellets, etc.

The NIR system provided to clients is based on the Bruker Matrix FT-NIR with an integrating sphere and PbS detector. The Bruker Matrix system is designed to work in a production environment, as it is waterproof, dust-proof and very rugged. An industrial notebook computer with touch screen function is also integrated with the NIR system. Via the Internet, it is possible to remotely monitor the instrument performance, and often to identify and solve problems remotely.

The system's use of Fourier transform (FT) NIR technology allows calibration models from any type of NIR system to be transferred since FT-NIR can be adjusted to simulate the spectrum of any type of NIR system, including a different brand of FT-NIR, dispersive NIR, AOTF NIR and diode array NIR. The wavelength accuracy

of FT-NIR also makes the Internet-enabled NIR system easier to use the same calibration model for different client NIR units.

**CENTRALIZED MODEL BUILDING**

QTA also builds single models for entire networks of NIR units, not just calibration models for individual NIR analyzers. It is very important that these models predict consistent results of the same sample from different client NIR analyzers. QTA uses special calibration techniques to compensate for instrument variances of the FT-NIR systems, so the calibration model can be used by different NIR units including the units already in service.

QTA models not only compensate for hardware variances, but also for variances from sample and environmental changes. Therefore, the same sample measured by the NIR instruments at different locations, different days or after hardware changes (such as light source and laser replacement) still produce consistent results.

Since all the models are stored at the central processor, they can be remotely developed or modified by NIR experts at the central location and used by all the client NIR systems simultaneously. Therefore, even without on-site NIR experts, the NIR applications for a client system can still be easily expanded.

Rapid, accurate, reliable on-site NIR analyses of multiple traits or material characteristics of both solids and liquids — without the need for skilled analytical personnel and without significant site-to-site, instrument-to-instrument, or day-to-day variations — is now a reality. And since all the spectra and predicted results are stored in a database that resides at the central location, either the current analyzed data or the historic data can be shared across the entire supply chain. **WG**

E-archive #12345. Find this story in the Article Archives at [www.World-Grain.com](http://www.World-Grain.com). Send comments and inquiries about this article to [worldgrain@sosland.com](mailto:worldgrain@sosland.com).

**QTA step-by-step**

1. Select analysis

2. Enter sample information

3. View results

**Superior Grain Management**

From temp cables and handheld monitoring to automated control systems, OPI offers the most advanced grain management system on the market. Call us today to turn your storage into a profit center.

North America: **1-800-661-1055** • International: **403-219-3177**  
**For more information, email: [opisystems@opisystems.com](mailto:opisystems@opisystems.com)**

**1&O StorMax™**  
[www.opistormax.com](http://www.opistormax.com)